

## Immunofluorescence on tobacco BY-2 suspension cells (GA fix)

cells should be happy and healthy.

*total time needed*

take approx. 700  $\mu$ l cell culture, add 10  $\mu$ l glutaraldehyde (70%) **in the hood**

incubate for 15 min on “labquake” rotator.

0.25 h

wash 3x with PBS (**in the hood**)

0.75 h

*how to wash:* spin 2 min @ 250 rpm w/ “soft stop” in Eppendorf microfuge ( $\approx 600 \times g$ ) – this may have to be repeated after the tubes have been turned 180° to get all cells down; carefully remove supernatant w/ pasteur pipette – remember to discard into GA waste; fill tube with fresh buffer; on rotator for  $\sim 5$  minutes

add washed cells to 0.1% NaBH<sub>4</sub> (sodium borohydride) in PBS and let sit over night in the hood

NaBH<sub>4</sub> has to be prepared fresh every time (best as the cells are being fixed). This step is supposed to remove the glutaraldehyde fluorescence

next morning:

PEI coat wells on microscope slide:

add 10  $\mu$ l of 0.1% polyethyleneimine per well, incubate in moist chamber for  $\geq 1$  h. (‘moist chamber’ is a petridish with wet filter paper – keep bottom of slide dry by propping it up on toothpicks, straws, etc) let water evaporate before loading cells (the draft at the chemical hood helps to speed this up)

remove supernatant and discard (sink), transfer 400  $\mu$ l cells w/ wide-bore tip to microfuge tube

have cells at  $\sim 3$ x packed volume in PBS, add 1:4 volume (i.e. 100  $\mu$ l) 5x enzymes, incubate for **15 minutes** (keep this volume  $\ll 500$   $\mu$ l so the following dilution/wash step is more efficient)

dilute by filling tube w/ PBS, wash 2x with PBS (do ‘double spin’ to ensure individual cells aren’t lost in the supernatant)

1h

let cells settle onto PEI-coated multi-well slide (use 12 $\mu$ l/well for 12-well slide) for 20 to 30 minutes

wash cells 1x with PBS

*how to wash on the slide:* carefully remove  $\sim 10$   $\mu$ l (i.e. less than you loaded) of the supernatant from the cells – pipette slowly, making sure that most cells stay put and **cells are always covered with buffer!** Don’t ever let them dry out! If you see little humps on the buffer surface they are already starting to plasmolyze! Add another 10  $\mu$ l buffer (PBS) to the cells and incubate for 5 minutes.

add  $\sim 10$   $\mu$ l 1% Nonidet (in PBS) for 15 minutes

wash 3x with PBS

2 h

add primary antibody (in PBS) for 0.5 to 1 h [for MTs use “4A1” anti-tubulin 1:40]

3h

wash 3x with PBS

3.5 h

add secondary antibody (in PBS) for 1 h [for MTs use Alexa-a-mouse at 1:400]

4.5 h

wash (at least) 3x with PBS

5 h

remove most of the supernatant (but leave cells covered!), add 1-2  $\mu$ l PBS/glycerol (1:1), view [fluorescence is stable for  $\geq 1$  week if coverslip is sealed and slide kept in dark (RT)]

## **Buffers etc:**

### **PBS (Phosphate Buffered Saline):**

(see Maniatis)

Keep it sterile so no bugs will grow.

### **5x Enzyme Mix**

use cellulase (Worthington) and pectinase (Fluka) at a final concentration of 1%  
to prepare 5x mix weigh out ~10 mg each of the enzymes and add ~200  $\mu$ l MSBG  
always prepare enzymes fresh (usually while washing PFA out)